



## FlowLine Bipore<sup>®</sup>

### ePTFE Vascular Grafts

- Unique Bipore design with two different fibril lengths for low thrombogenicity
- An additional ePTFE wrap guarantees enhanced suture retention and high burst strength
- Unique guideline indicates diameter and wall thickness of the graft
- Helical reinforcement for high resistance against kinking and compression
- Excellent and pliable handling and suture behaviour
- Simple and easy removal of the spiral reinforcement

# FlowLine Bipore®

## ePTFE Vascular Grafts

### Ordering information

Quick and reliable delivery in Europe within 24 hours.

Call our central office: +49 (0) 74 71 / 922 - 0

#### Thin Wall

#### Thin Wall with spiral reinforcement

Length (cm)	Diameter (mm)	Catalogue No.	Catalogue No.
10	5	10TW1005N	10TW1005S
20	6	10TW2006N	-
20	7	10TW2007N	-
20	8	10TW2008N	-
40	5	10TW4005N	10TW4005S
40	6	10TW4006N	-
50	6	10TW5006N	10TW5006S
50	7	10TW5007N	10TW5007S
50	8	10TW5008N	10TW5008S
70	5	10TW7005N	10TW7005S
80	6	10TW8006N	10TW8006S
80	7	10TW8007N	10TW8007S
80	8	10TW8008N	10TW8008S

#### Standard Wall

#### Standard Wall with spiral reinforcement

Length (cm)	Diameter (mm)	Catalogue No.	Catalogue No.
20	6	10SW2006N	-
20	7	10SW2007N	-
40	5	10SW4005N	-
50	6	10SW5006N	10SW5006S
50	7	10SW5007N	10SW5007S
50	8	10SW5008N	10SW5008S
70	5	10SW7005N	10SW7005S
80	6	10SW8006N	10SW8006S
80	7	10SW8007N	10SW8007S
80	8	10SW8008N	10SW8008S

# JOTEC®

SOLUTIONS FOR VASCULAR DISEASE

#### JOTEC GmbH

Lotzenäcker 23

D-72379 Hechingen

Tel. +49 (0) 74 71 / 922 - 0

Fax +49 (0) 74 71 / 922 - 100

www.jotec.net

#### JOTEC AG

Luzernerstr. 91

CH-5630 Muri

Tel. +41 (0) 56 670 91 - 00

Fax +41 (0) 56 670 91 - 02

